

CAPSULES	DESCRIPTION	TECHNOLOGY HEAD
Hyaluronic	Contains Cristalhyal that acts as a “molecular sponge” and thereby increases skins capacity to hold water – keeping it moist, nourished and healthy. It increases overall ammunity, moisturizes and prevents dehydration and repairs the skin by stimulation collagen and elastin synthesis	COMBINATION SKIN - <b>RED LIGHT</b> NORMAL SKIN - MICROCURRENT
BTX	Argireline® peptide solution C works by modulating muscle contraction to inhibit the production of wrinkles. It reduces wrinkles, smoothes skin and condition skin.	<b>MICROCURRENT</b>
Collagen	Mari Coll, derived from deep water fish, is a highly pure collagen that works to reinforce the skins immune system and fiht free radicals to slow down skin ageing and cell degraation. It moisturizes and nourishes the skin, protects the skin from harmful radiation and slows down skin ageing and cell degradation.	AGE PREVENTION - MICROCURRENT DRY SKIN – <b>RED LIGHT</b>
Natural Line Anti-Ageing	Pure Jojoba Oil natural sebum contains wax esters that act as a natural moisturizer and barrier for the skin. It is used for nursing skin and prevent dehydration.	NORMAL SKIN - <b>RED LIGHT</b> AGE PREVENTION - MICROCURRENT
Natural Rejuvenation	Shea Butter penetrates deep into the skin. It rejuvenates damaged cells, Restores elasticity and skin tone and protects and regenerates the skin.	NORMAL / COMBINATION / DRY - <b>RED LIGHT</b> SENSITIVE SKIN – <b>BLUE LIGHT</b>
Natural Line Peeling Exfoliation	A natural scrub, Sweet Almond Oil enriched with cellulose peeling with vitamins E & A. Removes Dead Skin Cells, Increase blood circulation and encourage cell turnover.	
Natural Line Calming Milk	Calendula Oil has antiseptic, anti fungal, anti-inflammatory and antibacterial properties. It soothes the skin, treats dry or chapped skin and treats burns, skin rashes, inflammation and wounds.	<b>BLUE LIGHT</b>
Active Prolong Whitening	Fuses 2 active whitening ingredients to lighten and even skin tone. It blocks epidermal melanin biosynthesis, minimises liver spots and even and lighening skin tone	AGE 16 – 35 - <b>BLUE LIGHT</b> AGE 36 – 46+ - MICROCURRENT
Active Prolong Repairing	Based on 2 different kinds of Hyaluronic Cirstalhyal and Primalhyal acids which together penetrate and fill skin tissues. It Prevents skin Hydration and moisturises, improves skin elasticity and firmness, reduction of wrinkles and provides cellular repair and wound healing	NORMAL / COMBINATION SKIN – MICROCURRENT  SENSITIVE SKIN – <b>BLUE LIGHT</b>
Active Prolong Neo Energy	Based on Neodermyl®: “needle-free” collagen and elastin filler provides a pure source of bio-energy for ski cells. It moisturises and nourishes the skin, improves skin elasticity and skin firmness and reduction of deep wrinkles.	<b>MICROCURRENT</b>